





ACTR3B Antibody

Product Code	CSB-PA878955ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1U1
Immunogen	Recombinant Human Actin-related protein 3B protein (189-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 3B (ARP3-beta) (Actin-like protein 3B) (Actin-related protein ARP4), ACTR3B, ARP11 ARP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTR3B


Image



Western blot
All lanes: ACTR3B antibody at 8µg/ml + HepG2 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 48, 38, 40 kDa
Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA878955ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA878955ESR1HU at dilution of 1:100