





ACTR3B Antibody

Product Code	CSB-PA878955ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1U1
Immunogen	Recombinant Human Actin-related protein 3B protein (189-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 3B (ARP3-beta) (Actin-like protein 3B) (Actin-related protein ARP4), ACTR3B, ARP11 ARP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTR3B


Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA878955ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA878955ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA878955ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.