




CD3D Antibody

Product Code	CSB-PA004930ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04234
Immunogen	Recombinant Human T-cell surface glycoprotein CD3 delta chain protein (22-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD3 delta chain (T-cell receptor T3 delta chain) (CD antigen CD3d), CD3D, T3D
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD3D

Image


Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004930ESR1HU at dilution of 1: 100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.