




CYGB Antibody

Product Code	CSB-PA006376ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWM9
Immunogen	Recombinant Human Cytoglobin protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoglobin (Histoglobin) (HGb) (Stellate cell activation-associated protein), CYGB, STAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYGB

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA006376ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.