




DHCR24 Antibody

Product Code	CSB-PA618005ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15392
Immunogen	Recombinant Human Delta(24)-sterol reductase protein (227-516AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delta(24)-sterol reductase (EC 1.3.1.72) (24-dehydrocholesterol reductase) (3-beta-hydroxysterol delta-24-reductase) (Diminuto/dwarf1 homolog) (Seladin-1), DHCR24, KIAA0018
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DHCR24

Image



Western blot


All lanes: DHCR24 antibody at 2µg/ml + Jurkat whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61, 55 kDa

Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA618005ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA618005ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.