




DIRAS3 Antibody

Product Code	CSB-PA006907ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95661
Immunogen	Recombinant Human GTP-binding protein Di-Ras3 protein (1-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein Di-Ras3 (Distinct subgroup of the Ras family member 3) (Rho-related GTP-binding protein Rho1), DIRAS3, ARHI NOEY2 RHOI
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	DIRAS3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA006907ESR2HU at dilution of 1:100