



MOAP1 Antibody

Product Code	CSB-PA014693ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BY2
Immunogen	Recombinant Human Modulator of apoptosis 1 protein (1-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modulator of apoptosis 1 (MAP-1) (MAP1) (Paraneoplastic antigen Ma4), MOAP1, PNMA4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MOAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.