



SIRT6 Antibody

Product Code	CSB-PA854057ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6T7
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-6 protein (1-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-6 (EC 3.5.1.-) (Regulatory protein SIR2 homolog 6) (SIR2-like protein 6), SIRT6, SIR2L6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIRT6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.