



SIRT6 Antibody

Product Code	CSB-PA854057ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6T7
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-6 protein (1-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-6 (EC 3.5.1.-) (Regulatory protein SIR2 homolog 6) (SIR2-like protein 6), SIRT6, SIR2L6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIRT6