





BCAT2 Antibody

Product Code	CSB-PA002601ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15382
Immunogen	Recombinant Human Branched-chain-amino-acid aminotransferase, mitochondrial protein (20-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Branched-chain-amino-acid aminotransferase, mitochondrial (BCAT(m)) (EC 2.6.1.42) (Placental protein 18) (PP18), BCAT2, BCATM BCT2 ECA40
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BCAT2


Image



Western blot
 All lanes: BCAT2 antibody at 14µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 45, 34 kDa
 Observed band size: 34 kDa



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA002601ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA002601ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.