




CD33 Antibody

Product Code	CSB-PA004925ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20138
Immunogen	Recombinant Human Myeloid cell surface antigen CD33 protein (49-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloid cell surface antigen CD33 (Sialic acid-binding Ig-like lectin 3) (Siglec-3) (gp67) (CD antigen CD33), CD33, SIGLEC3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD33

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA004925ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.