



TXNDC12 Antibody

Product Code	CSB-PA025371HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95881
Immunogen	Recombinant Human Thioredoxin domain-containing protein 12 protein (27-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin domain-containing protein 12 (EC 1.8.4.2) (Endoplasmic reticulum resident protein 18) (ER protein 18) (ERp18) (Endoplasmic reticulum resident protein 19) (ER protein 19) (ERp19) (Thioredoxin-like protein p19) (hTLP19), TXNDC12, TLP19
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TXNDC12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.