





CYP2F1 Antibody

Product Code	CSB-PA006426ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24903
Immunogen	Recombinant Human Cytochrome P450 2F1 protein (100-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 2F1 (EC 1.14.14.1) (CYP11F1), CYP2F1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CYP2F1

Image



Western blot
 All lanes: Cytochrome P450 2F1 antibody at 8µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 56, 37 kDa
 Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA006426ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.