



## stxB Antibody

<b>Product Code</b>	CSB-PA300809HA01EKZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P69179
<b>Immunogen</b>	Recombinant Enterobacteria phage H19B Shiga-like toxin 1 subunit B protein (21-89AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Enterobacteria phage H19B
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Shiga-like toxin 1 subunit B (SLT-1 B subunit) (SLT-1b) (SLT-1b) (Verocytotoxin 1 subunit B) (Verotoxin 1 subunit B), stxB, sltB
<b>Immunogen Species</b>	Enterobacteria phage H19B (Bacteriophage H19B)
<b>Research Area</b>	Others
<b>Target Names</b>	stxB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.