




ANXA4 Antibody

Product Code	CSB-PA001845ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09525
Immunogen	Recombinant Human Annexin A4 protein (110-319AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A4 (35-beta calcimedlin) (Annexin IV) (Annexin-4) (Carbohydrate-binding protein p33/p41) (Chromobindin-4) (Endonexin I) (Lipocortin IV) (P32.5) (PP4-X) (Placental anticoagulant protein II) (PAP-II) (Protein II), ANXA4, ANX4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANXA4

Image



Western blot

All lanes: Annexin A4 antibody at 4µg/ml

Lane 1: 293T whole cell lysate

Lane 2: HeLa whole cell lysate


Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 36, 28 kDa

Observed band size: 36 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA001845ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.