




RRM2B Antibody

Product Code	CSB-PA745339ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7LG56
Immunogen	Recombinant Human Ribonucleoside-diphosphate reductase subunit M2 B protein (1-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonucleoside-diphosphate reductase subunit M2 B (EC 1.17.4.1) (TP53-inducible ribonucleotide reductase M2 B) (p53-inducible ribonucleotide reductase small subunit 2-like protein) (p53R2), RRM2B, P53R2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RRM2B

Image



Western blot

All lanes: RRM2B antibody at 3µg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: HeLa whole cell lysate


Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 41, 35, 16, 8, 5, 49 kDa

Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA745339ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.