




HOXB1 Antibody

Product Code	CSB-PA010660ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14653
Immunogen	Recombinant Human Homeobox protein Hox-B1 protein (35-180AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-B1 (Homeobox protein Hox-2I), HOXB1, HOX2I
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXB1

Image



Western blot
 All lanes: Homeobox protein Hox-B1 antibody at 9µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 33, 25 kDa
 Observed band size: 33 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.