




PTPN1 Antibody

Product Code	CSB-PA019024DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18031
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 1 protein (281-435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 1 (EC 3.1.3.48) (Protein-tyrosine phosphatase 1B) (PTP-1B), PTPN1, PTP1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN1

Image



"Western blot
 All lanes: PTPN1 antibody at 2µg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: A431 whole cell lysate
 Lane 3: PC-3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 50 kDa
 Observed band size: 50 kDa
 "

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.