





TCOF1 Antibody

Product Code	CSB-PA619765DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13428
Immunogen	Recombinant Human Treacle protein (1189-1488AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Treacle protein (Treacher Collins syndrome protein), TCOF1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TCOF1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA619765DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA619765DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.