





COQ8A Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA850918ESR1HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NI60 |
| Immunogen | Recombinant Human Atypical kinase ADCK3, mitochondrial protein (1-220AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Atypical kinase COQ8A, mitochondrial (EC 2.7.-.-) (Chaperone activity of bc1 complex-like) (Chaperone-ABC1-like) (Coenzyme Q protein 8A) (aarF domain-containing protein kinase 3), COQ8A, ADCK3 CAB31 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | COQ8A |

Image



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA850918ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA850918ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.