



Pectate lyase 1 Antibody

Product Code	CSB-PA305864HA01JAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81294
Immunogen	Recombinant Juniperus ashei Pectate lyase 1 protein (22-367AA)
Raised In	Rabbit
Species Reactivity	Juniperus ashei
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Jun a 1) (allergen Jun a 1)
Immunogen Species	Juniperus ashei
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.