




SERPINB5 Antibody

Product Code	CSB-PA021073ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36952
Immunogen	Recombinant Human Serpin B5 protein (21-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin B5 (Maspin) (Peptidase inhibitor 5) (PI-5), SERPINB5, PI5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SERPINB5

Image



Western blot

All lanes: Serpin B5 antibody at 10µg/ml

Lane 1: HepG2 whole cell lysate

Lane 2: HeLa whole cell lysate


Lane 3: A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43, 26 kDa

Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded

human breast cancer using CSB-

PA021073ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.