



toxB Antibody

Product Code	CSB-PA321589HA01DUN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18177
Immunogen	Recombinant Peptoclostridium difficile Toxin B protein (1751-2366AA)
Raised In	Rabbit
Species Reactivity	Peptoclostridium difficile
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toxin B (EC 3.4.22.-), toxB, tcdB
Immunogen Species	Peptoclostridium difficile
Research Area	Others
Target Names	toxB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.