




PPIL4 Antibody

Product Code	CSB-PA851964LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUA2
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase-like 4 protein (1-492AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase-like 4 (PPIase) (EC 5.2.1.8) (Cyclophilin-like protein PPIL4) (Rotamase PPIL4), PPIL4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPIL4

Image



Western blot

All lanes: PPIL4 antibody at 4µg/ml

Lane 1: Hela whole cell lysate

Lane 2: Mouse kidney tissue


Lane 3: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58 kDa

Observed band size: 70 kDa



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA851964LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.