



tar Antibody, HRP conjugated

Product Code	CSB-PA356999HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07017
Immunogen	Recombinant Escherichia coli Methyl-accepting chemotaxis protein II protein (212-553AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methyl-accepting chemotaxis protein II (MCP-II) (Aspartate chemoreceptor protein), tar, cheM
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	tar
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.