



SIRPA Antibody, Biotin conjugated

Product Code	CSB-PA021334LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78324
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type substrate 1 protein (31-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type substrate 1 (SHP substrate 1) (SHPS-1) (Brain Ig-like molecule with tyrosine-based activation motifs) (Bit) (CD172 antigen-like family member A) (Inhibitory receptor SHPS-1) (Macrophage fusion receptor) (MyD-1 antigen) (Signal-regulatory protein alpha-1) (Sirp-alpha-1) (Signal-regulatory protein alpha-2) (Sirp-alpha-2) (Signal-regulatory protein alpha-3) (Sirp-alpha-3) (p84) (CD antigen CD172a), SIRPA, BIT MFR MYD1 PTPNS1 SHPS1 SIRP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SIRPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.