




HSPB8 Antibody

Product Code	CSB-PA891531ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJY1
Immunogen	Recombinant Human Heat shock protein beta-8 protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-8 (HspB8) (Alpha-crystallin C chain) (E2-induced gene 1 protein) (Protein kinase H11) (Small stress protein-like protein HSP22), HSPB8, CRYAC E2IG1 HSP22
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB8

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA891531ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.