



P4HTM Antibody

Product Code	CSB-PA878913LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXG6
Immunogen	Recombinant Human Transmembrane prolyl 4-hydroxylase protein (82-563AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane prolyl 4-hydroxylase (P4H-TM) (EC 1.14.11.-) (Hypoxia-inducible factor prolyl hydroxylase 4) (HIF-PH4) (HIF-prolyl hydroxylase 4) (HPH-4), P4HTM, PH4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	P4HTM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.