



SSX2 Antibody, HRP conjugated

Product Code	CSB-PA613674LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16385
Immunogen	Recombinant Human Protein SSX2 protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein SSX2 (Cancer/testis antigen 5.2) (CT5.2) (Synovial sarcoma, X breakpoint 2) (Tumor antigen HOM-MEL-40), SSX2; SSX2B, SSX2A;
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SSX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.