



hup1 Antibody, Biotin conjugated

Product Code	CSB-PA358373DD01FQI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A3H6
Immunogen	Recombinant Streptomyces lividans DNA-binding protein HU 1 protein (1-93AA)
Raised In	Rabbit
Species Reactivity	Streptomyces lividans
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein HU 1 (HSI), hup1, hup
Immunogen Species	Streptomyces lividans
Research Area	Others
Target Names	hup1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.