



TCEANC2 Antibody, HRP conjugated

Product Code	CSB-PA842718LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MN5
Immunogen	Recombinant Human Transcription elongation factor A N-terminal and central domain-containing protein 2 protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription elongation factor A N-terminal and central domain-containing protein 2, TCEANC2, C1orf83
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCEANC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.