



ANGPTL2 Antibody, Biotin conjugated

Product Code	CSB-PA891539LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKU9
Immunogen	Recombinant Human Angiotensin-converting enzyme 2 protein (23-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiotensin-converting enzyme 2 (Angiotensin-like protein 2), ANGPTL2, ARP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANGPTL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.