



nrdA Antibody, Biotin conjugated

Product Code	CSB-PA355522HD01EVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00452
Immunogen	Recombinant Escherichia coli Ribonucleoside-diphosphate reductase 1 subunit alpha protein (1-761AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonucleoside-diphosphate reductase 1 subunit alpha (EC 1.17.4.1) (Protein B1) (Ribonucleoside-diphosphate reductase 1 R1 subunit) (Ribonucleotide reductase 1), nrdA, dnaF
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	nrdA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.