



nrdA Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA355522HD01EVN |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P00452 |
| Immunogen | Recombinant Escherichia coli Ribonucleoside-diphosphate reductase 1 subunit alpha protein (1-761AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ribonucleoside-diphosphate reductase 1 subunit alpha (EC 1.17.4.1) (Protein B1) (Ribonucleoside-diphosphate reductase 1 R1 subunit) (Ribonucleotide reductase 1), nrdA, dnaF |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | nrdA |