



UBD Antibody, FITC conjugated

Product Code	CSB-PA025435LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15205
Immunogen	Recombinant Human Ubiquitin D protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin D (Diubiquitin) (Ubiquitin-like protein FAT10), UBD, FAT10
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	UBD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.