



Major allergen Api g 1, isoallergen 2 Antibody, HRP conjugated

Product Code	CSB-PA310896LB01DNL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P92918
Immunogen	Recombinant Apium graveolens Major allergen Api g 1, isoallergen 2 protein (1-159AA)
Raised In	Rabbit
Species Reactivity	Apium graveolens
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major allergen Api g 1, isoallergen 2 (Allergen Api g 1.0201) (allergen Api g 1)
Immunogen Species	Apium graveolens
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.