



PTGDS Antibody, HRP conjugated

Product Code	CSB-PA018969HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41222
Immunogen	Recombinant Human Prostaglandin-H2 D-isomerase protein (23-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostaglandin-H2 D-isomerase (EC 5.3.99.2) (Beta-trace protein) (Cerebrin-28) (Glutathione-independent PGD synthase) (Lipocalin-type prostaglandin-D synthase) (Prostaglandin-D2 synthase) (PGD2 synthase) (PGDS) (PGDS2), PTGDS, PDS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTGDS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.