



# HABP2 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA12159B0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q14520
<b>Immunogen</b>	Recombinant Human Hyaluronan-binding protein 2 protein (314-560AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Hyaluronan-binding protein 2 (EC 3.4.21.-) (Factor VII-activating protease) (Factor seven-activating protease) (FSAP) (Hepatocyte growth factor activator-like protein) (Plasma hyaluronan-binding protein) [Cleaved into: Hyaluronan-binding protein 2 50 kDa heavy chain; Hyaluronan-binding protein 2 50 kDa heavy chain alternate form; Hyaluronan-binding protein 2 27 kDa light chain; Hyaluronan-binding protein 2 27 kDa light chain alternate form], HABP2, HGFAL PHBP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	HABP2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.