



MCTS1 Antibody, HRP conjugated

Product Code	CSB-PA892143HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULC4
Immunogen	Recombinant Human Malignant T-cell-amplified sequence 1 protein (1-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malignant T-cell-amplified sequence 1 (MCT-1) (Multiple copies T-cell malignancies), MCTS1, MCT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MCTS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.