



SPT16 Antibody, FITC conjugated

Product Code	CSB-PA022955LC01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O82491
Immunogen	Recombinant Arabidopsis thaliana FACT complex subunit SPT16 protein (890-1074AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FACT complex subunit SPT16 (Facilitates chromatin transcription complex subunit SPT16), SPT16
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	SPT16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.