



fabD Antibody, HRP conjugated

Product Code	CSB-PA686101LB01FLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HGK3
Immunogen	Recombinant Staphylococcus aureus Malonyl CoA-acyl carrier protein transacylase protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malonyl CoA-acyl carrier protein transacylase (MCT) (EC 2.3.1.39), fabD
Immunogen Species	Staphylococcus aureus (strain COL)
Research Area	Others
Target Names	fabD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.