



fabD Antibody, Biotin conjugated

Product Code	CSB-PA686101LD01FLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HGK3
Immunogen	Recombinant Staphylococcus aureus Malonyl CoA-acyl carrier protein transacylase protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malonyl CoA-acyl carrier protein transacylase (MCT) (EC 2.3.1.39), fabD
Immunogen Species	Staphylococcus aureus (strain COL)
Research Area	Others
Target Names	fabD