



Casp12 Antibody, HRP conjugated

Product Code	CSB-PA004545LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O08736
Immunogen	Recombinant Mouse Caspase-12 protein (1-92AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caspase-12 (CASP-12) (EC 3.4.22.-), Casp12
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Casp12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.