



PTGES3 Antibody

Product Code	CSB-PA613587ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15185
Immunogen	Recombinant Human Prostaglandin E synthase 3 protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostaglandin E synthase 3 (EC 5.3.99.3) (Cytosolic prostaglandin E2 synthase) (cPGES) (Hsp90 co-chaperone) (Progesterone receptor complex p23) (Telomerase-binding protein p23), PTGES3, P23 TEBP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PTGES3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.