





APLP2 Antibody

Product Code	CSB-PA001912DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06481
Immunogen	Recombinant Human Amyloid-like protein 2 protein (340-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amyloid-like protein 2 (APLP-2) (APPH) (Amyloid protein homolog) (CDEI box-binding protein) (CDEBP), APLP2, APPL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	APLP2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA001912DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA001912DSR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.