




# PRKG1 Antibody

<b>Product Code</b>	CSB-PA619797ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q13976
<b>Immunogen</b>	Recombinant Human cGMP-dependent protein kinase 1 protein (1-300AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	cGMP-dependent protein kinase 1 (cGK 1) (cGK1) (EC 2.7.11.12) (cGMP-dependent protein kinase I) (cGKI), PRKG1, PRKG1B PRKGR1A PRKGR1B
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	PRKG1

## Image



### Western blot

All lanes: cGMP-dependent protein kinase 1 antibody at 3μg/ml + Mouse lung tissue

### Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 77, 78, 45 kDa

Observed band size: 77 kDa