




PRKG1 Antibody

Product Code	CSB-PA619797ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13976
Immunogen	Recombinant Human cGMP-dependent protein kinase 1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	cGMP-dependent protein kinase 1 (cGK 1) (cGK1) (EC 2.7.11.12) (cGMP-dependent protein kinase I) (cGKI), PRKG1, PRKG1B PRKGR1A PRKGR1B
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PRKG1

Image



Western blot
 All lanes: cGMP-dependent protein kinase 1 antibody at 3µg/ml + Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 77, 78, 45 kDa
 Observed band size: 77 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.