




CHRNA5 Antibody

Product Code	CSB-PA005391ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30532
Immunogen	Recombinant Human Neuronal acetylcholine receptor subunit alpha-5 protein (23-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal acetylcholine receptor subunit alpha-5, CHRNA5, NACHRA5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRNA5

Image



Western blot
 All lanes: Neuronal acetylcholine receptor subunit alpha-5 antibody at 12µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 53 kDa
 Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.