




AIPL1 Antibody

Product Code	CSB-PA882159ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN9
Immunogen	Recombinant Human Aryl-hydrocarbon-interacting protein-like 1 protein (1-384AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aryl-hydrocarbon-interacting protein-like 1, AIPL1, AIPL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AIPL1

Image



Western blot


All lanes: Aryl-hydrocarbon-interacting protein-like 1 antibody at 5µg/ml + Rat brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 44, 37, 42, 41 kDa

Observed band size: 44 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA882159ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.