




BAP1 Antibody

Product Code	CSB-PA849763ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92560
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase BAP1 protein (460-729AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase BAP1 (EC 3.4.19.12) (BRCA1-associated protein 1) (Cerebral protein 6), BAP1, KIAA0272
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BAP1

Image



Western blot
All lanes: BAP1 antibody at 8µg/ml + Hela whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 80 kDa
Observed band size: 80 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.